

Kera Sit2Sit



- The Kera was designed in consultation with Physios and Occupational Therapists for comfort and safety.
- The Kera allows for easy transfers with a single caregiver. Quickly and easily transfer your loved one from their bed, chair or toilet to another seated position.
- The Kera is light and manoeuvrable with a small footprint and no electrical parts enabling it to be conveniently stored or transported in your home.
- Kera 80 is suitable for transferring lighter people between seats and manual wheelchairs. Ideal for children.
- Unique chestpad and backstrap combination securely lifts by the upper body, no standing or sitting strength required!
- Accessory range to adapt to user ability and size, including children.
- Low model seat heights: 14.9 - 21.6 ' "
- High model seat heights: 21.6 - 25.6 ' "

Code	Display Name	Total Length	Total Width	Safe Working Load	Unit Weight	Leg Gait Range
UK-LSS398200	Kera Sit2Sit 80 - Low Model	37.4 - 48.4 "	23.8 "	12.6st	5.5st	16.3 - 43.3 "
UK-LSS398210	Kera Sit2Sit 80 - Tall Model	37.4 - 48.4 "	23.8 "	12.6st	5.5st	16.3 - 43.3 "
UK-LSS398220	Kera Sit2Sit 120 - Low Model	37.4 - 59.8 "	23.8 "	18.9st	5.5st	16.3 - 43.3 "
UK-LSS398230	Kera Sit2Sit 120 - Tall Model	37.4 - 59.8 "	23.8 "	18.9st	5.5st	16.3 - 43.3 "