

Aspire Vogue Carbon Fibre Seat Walker / Rollator



- **Super lightweight** - Carbon fibre frame.
- Now supplied with bonus cane holder (ref. image)
- **Luxury look and feel** – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- **Lightweight and ergonomic design** with small turning circle for excellent manoeuvrability.
- **“AutoLocate”** - Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate.
- **Locks when folded** – makes transportation and storage easier.
- **Advanced locking handbrakes** – for increased user safety when walking and sitting.
- **Unbreakable TPE super soft wheels** – offers high durability and a smooth ride for increased comfort.
- **Removable utility bag** – simple and stylish to remove for shopping. Can hold up to 0.78 st.
- *Unit weight includes accessories, bag and backstrap
- **Frame weight is minimal functional weight

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight
UK-WAF705350	Seat Walker - Aspire	27"	24"	27" - 36"	23.5 st	13.7 lb	Carbon	9"	18"	22"	22"	4.9 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight
	Vogue - Carbon Fibre						Fibre					