

Aspire Overbed Table - Thermoform Recessed Top



- Recessed Thermoform top to prevent spills
- Lightweight mobile table positions overbed or chair
- Can be used for eating, drawing or other activities
- Flat top ideal for hospital or home use
- Fully height adjustable
- "C" base shape fitted with four castors enables close positioning to the bedside or chair
- Powder coated mild steel frame construction
- Cream standard laminate top with vinyl wrap around the edge, or "thermoform" cream top which is heat resistant

Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Material
UK-BEA028000	Aspire Overbed Table - Thermoform Recessed Top - Cream	17"	31"	29" - 43"	Thermoform Top