

## Aspire Stride Seat Walker



- Side folding walking frame, ideal for storage and transport.
- Aluminium frame, lightweight and easy to manoeuvre.
- Ergonomic handles.
- Light touch lock brakes.
- Push button handle height adjustment.
- Soft padded seat for improved sitting comfort.
- Multi-terrain castors (Front 11' Back 9').
- Removable utility bag.

### Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Width	Seat To Floor Height	Seat To Footrest
UK-WAF750150	Aspire Stride Seat Walker - Red	32"	24"	31" - 38"	21 st	19.8 lb	Red	19"	20"	20"
UK-WAF750160	Aspire Stride Seat Walker - Purple	32"	24"	31" - 38"	21 st	19.8 lb	Purple	19"	20"	20"