

## Aspire 440mm Aluminium Commode



- Ideal for in home and community use
- Lightweight Aluminium design
- Easy clean, padded seat, back and arms
- Tool free height adjustment allows the commode to be easily adapted to varying needs
- Swing back comfort arms assist transfers
- Bowl and carrier for improved infection control
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- 4x total locking, height adjustable stainless steel castors provide brake access from all angles.
- \*Self Propelled model also available

### Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

| Code         | Display Name                            | Safe Working Load | Unit Weight | Seat Width | Backrest Height |
|--------------|---|-------------------|-------------|------------|-----------------|
| UK-BTC068010 | Aluminium Commode - Attendant Propelled | 21.5 st           | 33.3 lb     | 17"        | 19"             |