

Aspire Vida Folding Wheelchair - Attendant Propelled



- The Aspire Vida is a super lightweight transport wheelchair ideal for quick trips to the shops, appointments and catching up with friends.
- Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support.
- The 21 st safe working load caters to a broad range of users.
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers.
- Contoured and padded seating improves comfort and postural support.
- Attendant operated hand brakes with parking lock.
- Simple and intuitive folding for easy transportation.
- Removable leg rests and folding back for compact storage.
- Rear storage pocket for essentials.
- Available as attendant propelled or self-propelled.
- Anterior Seat Height - 19'
- Posterior Seat Height - 18'
- Anti-tippers available (WCA848275)
- Lap belt available (WCA848300)

Warranty Information

- 5 Years - Frame 1 Year - Fabric & Components

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Frame Weight	Backrest Height
UK-MWS449874	Aspire Vida Folding Wheelchair -	39"	23.5"	36.5"	21 st	19.5 lb	Mustard	16.5"	16"	19.5 lb	17"

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Frame Weight	Backrest Height
	AP - Mustard - 16"										
UK-MWS449886	Aspire Vida Folding Wheelchair - AP - Olive - 20"	39"	27"	36.5"	21 st	20.5 lb	Olive	16.5"	20"	20.5 lb	17"