

## Verto 4HC Height Adjustable Shower Commode - Attendant Propelled



- Robust, lightweight stainless steel design
- Depth and height adjustable polyurethane backrest which is removable for cleaning purposes
- Epoxy powder coated frame
- Height adjustable swing-away legrests
- 5' tracking swivel castors
- Plastic pan/bowl runners
- Seat and bowl/pan carrier sold separately

### Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

| Code         | Display Name   | Safe Working Load | Seat Width | Backrest Height |
|--------------|--|-------------------|------------|-----------------|
| UK-BTC069490 | Verto Shower Commode - Swing Away Footrest - Attendant Propelled - 19" | 28 st             | 19"        | 19"             |
| UK-BTC069450 | Verto Shower Commode - Swing Away Footrest - Attendant Propelled - 22" | 28 st             | 22"        | -               |