

Verto Attendant propelled shower commode with weight bearing platform



- Adjustable seat height in 1' increments
- Depth and height adjustable polyurethane backrest which is removable for cleaning purposes
- Weight bearing footplate
- Robust, lightweight stainless steel construction
- Height adjustable swingaway legrests
- Height adjustable armrests
- *Seat sold separately

Warranty Information

- 5 Year

Code	Display Name	Total Width	Safe Working Load	Seat Width	Seat To Floor Height	Seat To Footrest
UK-BTC069495	Verto Shower Commode - Weight Bearing Platform - Attendant Propelled - 19"	22"	28 st	19"	20" - 25"	20" - 25"