

Aspire PostureFit Chair



PostureFit Configurable Lift Recline Chairs offer adjustable and adaptable seating and options designed to provide unparalleled postural support and comfort.

- **Head and Neck Support:** Whether you want to read, watch TV, or simply relax, the effortless electric headrest adjustment allows for convenient support and positioning of your head and neck which reduces strain.
- **Zero Gravity Tilt in space:** Creates a feeling of weightlessness and effectively redistributes sitting pressure away from vulnerable bony prominences, assists with oedema, and relaxes muscles.
- **Relieve Tired Legs and Enhance Circulation:** Independent legrest facilitates leg elevation with heel support. Recline to a comfortable position to reduce limb swelling, improve blood circulation, and alleviate fatigue.
- **Targeted Postural Support and Interchangeable Backrests:** All backrest cushions are zipper accessible and able to have filling or support blocks added or removed to customise the comfort and support to suit an individual's unique posture.
- **Air Cushion Technology:** Fitted within the seating support is a sectioned air cushion that provides targeted pressure relief and is segmented to ensure upright positioning is maintained.
- **After-Market Cushion Integration:** The layered seat structure accommodates after-market pressure care cushions such as ROHO replacing the air cushion where necessary, without compromising chair performance or infection control outcomes.
- **Medical Grade Foam:** High density foam allows for prolonged periods of sitting.
- **Integrated Footrest Extension:** Suits most user heights without the need to alter the legrest or height of the chair.
- **Tool-Free Seating Adjustment:** Therapists or carers can quickly and easily make adjustments without the need for any tools.
- **Effortless Mobility and Safety:** Directional castors are engaged by activating a rear

carer foot pedal. Castors are not in contact with the floor when in a seated position, ensuring stability and preventing accidental tipping or sliding.

- Transport mode: For safety, all electrical functions are inactivated while castors are engaged to prevent the chair being activated while mobile.

Warranty Information

- 2 Years Frame & Electronics
- 1 Year Handset & Cables

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
UK-CHP198250	Aspire PostureFit Chair - Small/Medium	30"-32"	44.9"	25st	18.7"-20.7"	19"-21"	19.3"	19.3"	29"
UK-CHP198260	Aspire PostureFit Chair - Tall Medium/Large	32.3"-34.3"	46.5"	25st	20.7"-22.6"	21"-23"	19.3"	19.3"	31"