



SOCIALITE WHEELCHAIR

ARRIVE IN STYLE

The Aspire Socialite is an ultra light weight transport chair perfect for short trips and outings with friends. The generous 115kg safe working load caters to a broad range of users.

Padded and breathable seat and back provide optimum comfort and postural support. Cut away arms improve table access and the offset uprights create hip space for coats and warmer clothing.

Less than
**9kg /
19.8lbs**



FEATURES & BENEFITS

- Sleek European design
- Stylish contemporary colours
- One of the lightest transport chairs available
- Perfect for outings and appointments
- Attendant operated hand brakes with parking lock
- Anti-tippers available – WCA848270
- Simple and intuitive folding for increased easy transportation
- Removable leg rests and folding back for compact storage
- Contoured seating improves comfort and postural support
- Rear storage pocket for essentials
- Available in attendant (12 in) and self-propelled (22 in) options



FOLDED MODE



Attendant Propelled



Self Propelled



UPHOLSTERY
Red/Black

FRAME
Red



UPHOLSTERY
Silver/Black

FRAME
Silver



SWL
115kg



SWL
18st

1 YEAR WARRANTY
Components

5 YEARS WARRANTY
Frame

SOCIALITE WHEELCHAIR

ATTENDANT PROPELLED
MWS449800 – SILVER
MWS449810 – RED

SELF PROPELLED
MWS449820 – SILVER
MWS449830 – RED

Overall Width	Seat Depth	Seat Width	Seat Height Anterior	Seat Height Posterior	Backrest Height	Armrest Height	Product Weight	Min. Chair Weight
mm/inches	mm/inches	mm/inches	mm/inches	mm/inches	mm/inches	mm/inches	kg/lbs	kg/lbs
595 / 23.4	430 / 17	450 / 17.7	475-500 / 18.7-19.7	450 / 17.7	410 / 16.1	200 / 7.9	9 / 19.8	8 / 17.6
675 / 26.6	430 / 17	450 / 17.7	475-500 / 18.7-19.7	450 / 17.7	410 / 16.1	200 / 7.9	13 / 28.7	8 / 17.6