



# ALUMINIUM SHOWER COMMODE



Self Propelled Model  
24" Lightweight Wheels

Attendant Propelled Model



View our Aspire Shower Commode Accessory Flyer

- **Height Adjustable, Swing Away Footrests** for optimal positioning and ease of transfers
- **4 x Total Locking, Height Adjustable Stainless Steel Castors** provide brake access from all angles
- **24" Lightweight Self propelling Wheels** with hand brakes and anti-tippers
- **Swing Back Comfort Arms** assist transfers
- **Tool-free Height Adjustment** allows the commode to be easily adapted to varying needs
- **Improved Infection Control** with bowl and carrier

Commode is registered in the ARTG AUST R ARTG 236829  
DVA contracted and AS/NZS approved

**5 YEAR WARRANTY**  
Corrosion

**1 YEAR WARRANTY**  
Components



Code	Description	Overall Width	Overall Height	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height
BTC068010	Attendant Propelled	460mm	990 – 1040mm	440mm	495mm	220mm	550mm	550mm
BTC068020	Self Propelled	610mm	970 – 995mm	440mm	490mm	235mm	525mm	510mm