## **ComfiMotion Collection**

Luxurious adjustable bed collection for an enhanced sleeping experience





We believe that creating a supportive and luxurious sleep environment can significantly improve health and wellbeing.

Our ComfiMotion Collection is designed to blend seamlessly into the home environment, whilst improving comfort and independence.


## Create Your ComfiMotion



**Bed Function** 



Headboard Style



**Mattress Support** 



Accessories

© Aidacare Pty Ltd 2024.


No part of this brochure may be copied, distributed, reproduced, scanned or stored on any electronic database, whether in whole or in part in any form or by any means

Due care has been taken to ensure all the information contained in this brochure is correct at the time of printing. All measurements, pictures, colours and weight capacities are to be used as a guide only. We reserve the right to modify the design or appearance of any product displayed in this brochure without prior notice.

ComfiMotion is a registered trademark of Aidacare Pty Ltd.

# **Activ Care**

### Ultra Low Adjustable Bed



### Supporting the modern, independent lifestyle.

The Activ Care bed combines the best of in-home care, functionality, positioning and low height safety with luxurious and contemporary style for enhanced comfort, care and wellbeing.

#### Features

- Ultra Low Height: Minimum platform height of 195mm reduces falls risk, improving user safety.
- Customisable Safe Entry/Exit Height: Programmable single button press ensures perfect bed height to maximise independence.
- Underbed Clearance: Provides clear access for patient lifters and overbed tables.
- Battery Backup: Provides confidence and maintains independence during power outages.
- Night Assist: Luminous underbed halo for safe entry and exit.
- Illuminated Handset: Large button interface is intuitive and highly visible at all times.
- Massage: Relaxing vibration modes for personalised neck, shoulder and back relief.

#### **Functions**

Hi-Lo Functionality: Assists with manual handling and helps to mitigate fall risks.

Head Tilt: Relaxes neck muscles, providing optimal support whilst reading or watching TV.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Auto-Regression in the Backrest: Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution.

Clinical Positioning: One touch raises back and legs to alleviate cardiac exertion, improving resting comfort.

Leg Raise: Promotes circulation, relieving tired, swollen, fatigued or painful limbs.

Memory Function: Allows the user to set and return to their preferred position at any time with one click of a button.

**Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.

Seated Position: Assists with bed entry and exit and upright dining.

Trendelenburg (and Reverse): Can assist care givers when repositioning a user\*.

Lockable Castors: Eight heavy duty dual castors.

#### ▲ Companion Sleeping

See Configurations section for more companion options.

**Upholstery** 



#### **Sizes** (Length x Width x Height Range)

**Long Single** 2030 x 920 x 195 - 700mm BEB046740

King Single 2030 x 1070 x 195 - 700mm BEB046750

Long Double 2030 x 1370 x 195 - 700mm BEB046760

Queen


2030 x 1520 x 195 - 700mm BEB046770

\*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.



# Care

### Adjustable Bed



#### For comfort you deserve.

The Care adjustable bed blends seamlessly into the home without compromising on the therapeutic needs of the user or safety of partners or carers. Select from multiple position options designed to maximise independent mobility and improve sleep quality.

#### **Features**

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- Illuminated Handset: Large button interface is intuitive and highly visible at all times.
- Enclosed Frame Design: Ensures optimal positioning when getting into bed and safely covers the moving base platform.
- Safe-Lowering Technology: Ensures that back, leg and Hi-Lo functions only lower under the force of gravity, reducing potential for accidents or injuries.
- Box-on-Box Design: Facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.

#### **Functions**

**Hi-Lo Functionality:** Assists with manual handling and helps to mitigate fall risks.

**Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

**Clinical Positioning:** Raises back and legs to alleviate cardiac exertion, improving resting comfort.

**Lockable Castors:** Four heavy duty castors.

Leg Raise: Promotes circulation, relieving tired, swollen, fatigued or painful limbs.

**Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.

**Seated Position:** Assists with bed entry and exit and upright dining.

**Trendelenburg (and Reverse):** Can assist care givers when repositioning a user\*.

### Variable Bed Height Options



Upholstery

\*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.



Low Height 280 – 580mm with Low Profile Bed Legs



Standard 350 – 650mm with Standard Castors



Accessibility 400 – 700mm with Large Castors (optional)



▲ Companion Sleeping

Sizes (Length x Width x Height Range)

#### **Long Single**

2030 x 920 x 280 - 580mm<sup>1</sup> BEB046720

#### **King Single**

2030 x 1070 x 280 – 580mm<sup>1</sup> BEB046730

1. With Low Profile Bed Legs.







SWL

Visit comfimotion.com.au for Limited Warranty T&Cs

# Relax

## Fixed Height Adjustable Bed



### Supporting the modern, independent lifestyle.

The Relax Bed has a fixed base and combines luxurious finishes with multiple positioning and therapeutic functions to support users and carers alike.

#### Features

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- Massage: Relaxing vibration modes for personalised neck, shoulder and back relief.
- Press & Hold Safety Handset: Intuitive, large button interface enhances usability.
- Box-on-Box Design: Facilitates easier entry and exit, allowing true flush positioning of beds in companion or dual adjustable bed configurations.
- German Engineered Motors: For smooth, quiet adjustment.
- Personalised Positioning: Hand control allows users to find their perfect position to maximise comfort and relaxation.
- USB Charge: Connect and charge your devices conveniently.

#### **Functions**

**Head Tilt:** Relaxes neck muscles, providing optimal support whilstreading or watching TV.

**Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

**Leg Raise:** Promotes circulation, relieving tired, swollen and/or painful limbs.

**Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.

**Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button.

#### **Upholstery**



#### **Sizes** (Length x Width x Height Range)

Long Single 2030 x 920 x 280 – 700mm BEB046900 King Single 2030 x 1070 x 280 – 700mm BEB046910 Long Double 2030 x 1370 x 280 – 700mm BEB046920 Long Queen 2030 x 1520 x 280 – 700mm BEB046930



▲ Companion Sleeping

See Configurations section for more companion options.



Visit comfimotion.com.au for Limited Warranty T&Cs

ID & O

# Mattresses

## Engineered for Motion

Our ComfiMotion Pocket Spring and Memory Foam mattresses are designed with the flexibility to follow the contouring profile of an adjustable bed.

### **Pocket Spring**

The ComfiMotion Pocket Spring mattress offers zoned support springs that respond and contour to the individual body shape to improve comfort, maximising support and sleep quality.



**Transitional Layered Pocket Springs** arranged in weight activated support zones.




Integrated Side-Walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.



Luxuriously Soft Pillowtop Layer provides a plush immersive feeling complimented by the responsive pocket springs.

YEAR LIMITED

WARRANTY



▲ Pocket Spring Mattress

## **Memory Foam - Classic**

(Length x Width x Height)

**Long Single Classic** 2030 x 920 x 200mm BEM469910

**King Single Classic** 

2030 x 1070 x 200mm BEM469920

**Long Double Classic** 

2030 x 1370 x 250mm BEM469922

2030 x 1520 x 250mm BEM469924

## Classic & Pressure Care

The ComfiMotion Classic and Pressure Care are memory foam mattresses that contour to the user's body, providing personalised comfort and transitional pressure relief. Integrated side-walls increase roll resistance at the bed edge and assist users, partners or carers with extra rigidity to ease bed entry and exit.



Contoured Immersion: Specially layered memory foam provides contoured immersion for support and pressure relief.



Superior Breathability: Up to 8x more breathable than conventional memory foams.



Anti-microbial: Environmentally safe treatment is added to the foam which reduces odours and kills 99.9% of bacteria.



**Pressure Care Incontinence Cover** 

Breathable, 4-way stretch waterproof cover to reduce shear, friction and maintain skin micro-climate.



Integrated Side–Walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.

### **Pocket Spring**

(Length x Width x Height)

SWL

180kg

LS & KS

**Pocket Spring Mattress – Long Single** 2030 x 920 x 240mm BEM469800

**Pocket Spring Mattress - King Single** 2030 x 1070 x 240mm BEM469810

**Pocket Spring Mattress - Long Double** 2030 x 1370 x 240mm BEM469820

Pocket Spring Mattress - Queen 2030 x 1520 x 240mm BEM469830

Queen Classic

#### **Memory Foam - Pressure Care** (Length x Width x Height)

**Long Single Pressure Care** 2030 x 920 x 200mm BEM469930

**King Single Pressure Care** 2030 x 1070 x 200mm BEM469940 **Long Double Pressure Care** 

2030 x 1370 x 250mm BEM469942 Queen Pressure Care

2030 x 1520 x 250mm BEM469944



PRESSURE CARE







11

Visit comfimotion.com.au for Limited Warranty T&Cs

360ka

Visit comfimotion.com.au for Limited Warranty T&Cs

### Clinical Considerations

**Accessibility and Inclusivity:** Unexposed bed surround with padded upholstery removes gaps between beds in co-sleeping arrangements and promotes safe entry and exit.

**User Independence:** Height adjustability improves user mobility and assists carers with manual handling.

**Flexible Positioning:** Caters to various user requirements such as respiratory conditions, snoring, sleep apnoea, pressure redistribution, joint or muscle pain, lymphoedema and swelling.

**Vision and Hand Function:** Handset with a clear and intuitive button display, assists users with reduced vision or hand function.

**Falls Prevention:** Integrated bed lighting and a minimum platform height of 195mm (Activ Care) and 280mm (Care and Relax) helps to reduce the likelihood and severity of falls in the home.

**User Safety:** Battery backup provides peace of mind in the event of a power outage (Standard on Activ Care, optional accessory on Care).

### Accessories



#### Crook Handle Bed Stick

Care/Relax Bed BEA830200 Activ Care Bed BEA830211 Provides multiple grip points to assist repositioning in bed, and safe entry and exit.



#### **Extendable Side Rail**

Function & Style Support

Activ Care Bed BEA830205

Provides a telescoping bed rail support for users at risk of falls.


### Headboards



#### Allure

Stylish freestanding headboard Subtle square padded design

Long Single BEA830120 King Single BEA830130 Long Double BEA830132 Queen BEA830134 King BEA830136



#### Chesterfield


Stylish freestanding headboard Chesterfield padded design

Long Single BEA830100
King Single BEA830110
Long Double BEA830112
Queen BEA830114
King BEA830116

#### **Bed Rail**

Care/Relax Bed BEA830210 Activ Care Bed BEA830212 Improves bed edge

Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility.



#### Split King Upgrade Kit

Activ Care Bed LS & KS only **BEA830230** 

Ensures simultaneous Hi–Lo platform movement whilst still allowing individual positioning requirements when beds are in partner configuration.

#### **Simplified Handset**

Activ Care Bed **BEA830214**For clients with reduced cognitive ability.



#### **Partner Handset**

Activ Care Bed BEA830213

Allows double and queen beds to be operated from either side of the bed.

### Pillows




#### **Breeze Cooling Gel Memory Foam**

- Adaptive and supportive sleep surface
- Helps regulate sleep surface temperature
- Accommodates side and back sleepers **BEP050100**



#### **Plush Memory Foam and Fibrefill**


- Familiar sleep surface
- Improved head and neck support
- Accommodates side and back sleepers **BEP050110**



### Large Accessibility Castors

Care Bed only (4 x required) **MWT451080** 

Larger castors raise the bed frame to 155mm, facilitating overbed table and patient lifter access.



#### **Battery Backup**

\*Standard on Activ Care Bed
BEA830220

Provides peace of mind in the event of a power outage.

### Collection Summary

Care Adjustable Bed – (LS)	BEB046720
Care Adjustable Bed – (KS)	BEB046730
Relax Adjustable Bed – (LS)	BEB046900
Relax Adjustable Bed – (KS)	BEB046910
Relax Adjustable Bed – (LD)	BEB046920

BEB046930

- 1
BEB046740
BEB046750
BEB046760
BEB046770
BEB046800

#### **Mattresses**

Relax Adjustable Bed - (Q)

Beds

Pocket Spring Mattress – (LS)	BEM469800
Pocket Spring Mattress – (KS)	BEM469810
Pocket Spring Mattress – (LD)	BEM469820
Pocket Spring Mattress – (Q)	BEM469830
Memory Foam Mattress – Classic (LS)	BEM469910
Memory Foam Mattress – Classic (KS)	BEM469920
Memory Foam Mattress – Classic (LD)	BEM469922
Memory Foam Mattress – Classic (Q)	BEM469924

Memory Foam Mattress – Pressure Care (LS) _ BEM469930
Memory Foam Mattress – Pressure Care (KS) _ BEM469940
Memory Foam Mattress – Pressure Care (LD) _ BEM469942
Memory Foam Mattress – Pressure Care (Q) BEM469944



### Headboards

Allure Headboard – (LS)	BEA830120
Allure Headboard – (KS)	BEA830130
Allure Headboard – (LD)	BEA830132
Allure Headboard – (Q)	BEA830134
Allure Headboard – (K)	BEA830136

Chesterfield Headboard – (LS)	BEA830100
Chesterfield Headboard – (KS)	BEA830110
Chesterfield Headboard – (LD)	BEA830112
Chesterfield Headboard – (Q)	BEA830114
Chesterfield Headboard – (K)	BEA830116



### **Pillows**

Breeze Pillow \_\_\_\_\_\_\_BEP05010

BEP050100 Plush Pillow



#### Accessories

Crook Handle Bed Stick – Care/Relax Bed	BEA830200
Crook Handle Bed Stick – Activ Care Bed	BEA830211
Bed Rail – Care/Relax Bed	BEA830210
Bed Rail – Activ Care Bed	BEA830212
Extendable Bed Rail	BEA830205

Large Accessibility Castors – Care Bed	MWT451080
Battery Backup – Care/Relax Bed	BEA830220
Split King Upgrade Kit – Activ Care Bed	BEA830230
Simplified Handset – Activ Care Bed	BEA830214
Partner Handset – Activ Care Bed	BEA830213

### Configurations for Companions

Whichever ComfiMotion package you select, rest assured there is a companion solution available to suit those wishing to sleep alongside any ComfiMotion Adjustable bed.

#### What to Order

**KING** 2030 x 1840



KING PLUS 2030 x 1990



**SUPERKING** 2030 x 2140



## Combine any two singles below:

Care Long Single BEB046720

Active Care Long Single

Fixed Long Single BEB046800
Relax Long Single BEB046900

### Combine any single with king single below:

Relax Long Single BEB046900
Care Long Single BEB046720

Active Care Long Single

Fixed Long Single

BEB046800

Relax King Single BEB046910
Care King Single BEB046730

Active Care King Single

## Combine any two king singles below:

Care King Single BEB046730

Active Care King Single

BEB046750

Relax King Single BEB046910

Long Double and Queen Beds can not be paired with other beds in the range.



#### **Dual Adjustable Beds**

Designed to sit flush when positioned together to enable close monitoring and enhanced care. This configuration offers the flexibility of independent adjustability, without compromising on the sense of inclusivity and intimacy between partners.

#### Split King Upgrade Kit

Ensures simultaneous Hi-Lo
ple platform movement whilst still
allowing individual positioning
requirements when beds are in
partner configuration.

**Split King Upgrade Kit** 

Compatible with active care bed

BEA830230

#### Fixed Companion Bed

Matching single bed without any adjustability designed to sit flush when positioned together. It facilitates a side-by-side sleeping arrangement at a fixed height of 370mm, without compromising on the users care and assistance requirements.

Fixed Companion Bed – (LS) 2030 x 920 x 370mm BEB046800

LS – Long Single

KS – King Single

LD - Long Double

**Q** – Queen

K – King





#### AIDACARE PTY LTD

Building 3A, 1 Moorebank Avenue, Moorebank NSW 2170 Australia

1300 133 120 | aidacare.com.au

#### AIDACARE LIMITED

Arcadia House, Maritime Walk, Ocean Village, Southampton SO14 3TL United Kingdom

01622 541235 | aidacare.co.uk

MADE IN PRC

ASBR28.96 040624

