



EVOKE 2 WHEELCHAIR

The Aspire EVOKE 2 and EVOKE 2 HD are highly adjustable, lightweight, aluminium wheelchairs designed to enable a personalised fit for every user. A comprehensive range of compatible accessories is available for increased independence and enhanced clinical outcomes.




CRASH TESTED
International Standard
ISO7176-19:2008

1 YEAR WARRANTY
Components

2 YEARS WARRANTY
Frame



MATT BLACK
FRAME



BLACK
UPHOLSTERY

EVOKE 2



EVOKE 2 HD



Swing-In / Out, Height Adjustable & Removable Foot Rests

This facilitates user transfers and allows optimal positioning for the changing needs of the user

Strong, Nylon Upholstery with a Comfortable Cushion & Padded Backrest

For increased comfort, a padded overlay is available as an accessory

Depth Adjustable Seating with Rear Angle Adjustability

For more accurate positioning, customisability and increased user comfort

Quick Release 24" Solid Rear Tyres

The wheels are easily removed to reduce the weight of the chair for transport and storage



ADJUSTMENTS



Seat & Backrest Height



Seat Depth



Backrest Angle

OPTIONAL ACCESSORIES



Contoured Padded Overlay*



Attendant Brake Kit & Conversion Kit**



*Scan here for more information on the Contoured Padded Overlay



**Scan here for more information on the Evoke 2 Wheelchair Attendant Brake options

EVOKE 2 WHEELCHAIR

	Seat Width in	Seat Depth in	Seat Height in	Overall Width in	Backrest Height in	Armrest Height in	Min. Chair Weight lbs	SWL st
MWS449505 Standard 14"	14	12 – 18	17 – 20	22	17 – 19	8 – 9	33	24
MWS449515 Standard 16"	16	12 – 18	17 – 20	24	17 – 19	8 – 9	33	24
MWS449525 Standard 18"	18	12 – 18	17 – 20	26	17 – 19	8 – 9	33	24
MWS449535 Standard 20"	20	16 – 22	17 – 20	28	17 – 19	8 – 9	33	24
MWS449545 HD 22"	22	16 – 22	17 – 20	31	17 – 19	8 – 9	35	30
MWS449555 HD 24"	24	16 – 22	17 – 20	32	17 – 19	8 – 9	36	30

ASFY127.6.210125