



# EVOKE 2 WHEELCHAIR

The Aspire EVOKE 2 and EVOKE 2 HD are highly adjustable, lightweight, aluminium wheelchairs designed to enable a personalised fit for every user. A comprehensive range of compatible accessories is available for increased independence and enhanced clinical outcomes.



CRASH TEST  
**APPROVED**  
AS/NZS 3696.19

**1** YEAR WARRANTY  
Components

**2** YEARS WARRANTY  
Frame



**MATT BLACK**  
FRAME



**BLACK**  
UPHOLSTERY

**EVOKE 2**



**EVOKE 2 HD**



**Swing-In / Out, Height Adjustable & Removable Foot Rests**

This facilitates user transfers and allows optimal positioning for the changing needs of the user

**Strong, Nylon Upholstery with a Comfortable Cushion & Padded Backrest**

For increased comfort, a padded overlay is available as an accessory

**Depth Adjustable Seating with Rear Angle Adjustability**

For more accurate positioning, customisability and increased user comfort

**Quick Release 24" Solid Rear Tyres**

The wheels are easily removed to reduce the weight of the chair for transport and storage

Pictured with optional transit conversion kit



\*Scan here for more information on the Evoke 2 Wheelchair Attendant Brake options

**ADJUSTMENTS**



**Seat & Backrest Height**



**Seat Depth**



**Backrest Angle**

**OPTIONAL ACCESSORIES**



**Contoured Padded Overlay**



**Attendant Brake Kit & Conversion Kit\***

**EVOKE 2 WHEELCHAIR**

	Seat Width in	Seat Depth in	Seat Height in	Overall Width in	Backrest Height in	Armrest Height in	Min. Chair Weight lbs	SWL st
MWS449505 Standard 14"	14	12 – 18	17 – 20	22	17 – 19	8 – 9	33	24
MWS449515 Standard 16"	16	12 – 18	17 – 20	24	17 – 19	8 – 9	33	24
MWS449525 Standard 18"	18	12 – 18	17 – 20	26	17 – 19	8 – 9	33	24
MWS449535 Standard 20"	20	16 – 22	17 – 20	28	17 – 19	8 – 9	33	24
MWS449545 HD 22"	22	16 – 22	17 – 20	31	17 – 19	8 – 9	35	30
MWS449555 HD 24"	24	16 – 22	17 – 20	32	17 – 19	8 – 9	36	30

ASFY127.3.011024