



VIDA WHEELCHAIR

The Aspire VIDA is an extremely lightweight wheelchair that can be lifted in and out of a vehicle easily by carers. The VIDA range comes in various fun, bold colours and is ideal for quick trips to the shops, appointments, and catching up with friends. The ergonomic backrest and seat shape mimic body contours for comfort and postural support. The high safe working load (21st) caters to a wide range of users.



UPHOLSTERY
Mustard/Black
FRAME
Grey
SEAT WIDTH
16"



UPHOLSTERY
Green/Black
FRAME
Silver
SEAT WIDTH
18"



UPHOLSTERY
Orange/Black
FRAME
Silver
SEAT WIDTH
18"



UPHOLSTERY
Olive/Black
FRAME
Grey
SEAT WIDTH
20"

Minimum
Chair Weight
from
7.5kg

SWL
21st

5 YEARS WARRANTY
Frame

1 YEAR WARRANTY
Components

FEATURES & BENEFITS

- Available as attendant propelled or self-propelled
- Three seat width sizes: 16", 18", and 20", each with its own unique colour
- The 16" seat width comes in Mustard
- The 18" seat width is available in Green and Orange
- The 20" seat width is available in Olive
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers
- Contoured and padded seating improves comfort and postural support
- Attendant operated hand brakes with parking lock
- Simple and intuitive folding for easy transportation
- Removable leg rests and folding back for compact storage
- Rear storage pocket for essentials
- Anti-tippers available (WCA848275)
- Lap belt available (WCA848300)



Attendant Propelled (AP) & Self Propelled (SP) Vida Specifications

16" Mustard



MWS449874 - AP
MWS449876 - SP

18" Green or Orange



MWS449855 - AP
MWS449860 - SP

20" Olive



MWS449886 - AP
MWS449888 - SP



VIDA	Overall Length in	Overall Width in	Overall Height in	Seat Depth in	Seat Width in	Anterior Seat Height in	Posterior Seat Height in	Backrest Height in	Armrest Height in	Min Chair Weight lbs	Rear Wheel in
Attendant Propelled											
MWS449874	39	23.5	36.5	16.5	16	19.5	18	17	8	19.5	13
MWS449855	39	25.5	36.5	16.5	18	19.5	18	17	8	19.5	13
MWS449850	39	25.5	36.5	16.5	18	19.5	18	17	8	19.5	13
MWS449886	39	27	36.5	16.5	20	19.5	18	17	8	20.5	13
Self Propelled											
MWS449876	41	25	37.5	16.5	16	19.5	18	17	8	16.5	22
MWS449860	41	27	37.5	16.5	18	19.5	18	17	8	16.5	22
MWS449865	41	27	37.5	16.5	18	19.5	18	17	8	16.5	22
MWS449888	41	28	37.5	16.5	20	19.5	18	17	8	17.6	22