

VOGUE ADVENTURE SEAT WALKER

Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super sized front wheels can manage any surface, perfect for the adventurous user!

COLOURS:



Gunmetal Grey



Wider Seat supports up to 31.5st



Large Multi-Terrain Wheels



Folds to 10"

1 YEAR WARRANTY

DANISH DESIGN Luxury Look & Feel Danish Designed

FEATURES & BENEFITS

- **Ergonomic Design** offers personalised height adjustment which promotes correct posture
- **Handbrakes** lock for user safety and stability
- **Simple** folding mechanism locks for easy transportation and storage
- **Compact** turning circle for improved mobility in tight spaces
- **Storage** removable bag carries up to 5kg of essentials
- **Bonus** Walking cane holder



ADVENTURE SEAT WALKER	Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size (Front/Rear)	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
GUNMETAL GREY WAF705455	30"	27.5"	35.5"-41"	9"	21"	21"	10"	11" / 9"	19.5lbs	22lbs	11lbs

** Includes accessories, bag and backstrap.