

VOGUE INDOOR WALKER

Ideal for use within the home, the Indoor Walker is perfect for portable dining and room to room mobility support.

COLOURS:



Red / Black



Violet / White



Folds to 11"

1 YEAR WARRANTY



FEATURES & BENEFITS

- **Adjustable Height:** Promotes correct posture by allowing personalised height adjustments.
- **Handbrakes for Stability:** Locks in place for user safety.
- **Folds Easily:** Simplifies transport and storage.
- **Safe Meal Transport:** Designed to carry meals and hot drinks securely.
- **Compact Turning:** Ideal for manoeuvring in tight spaces.
- **Removable Bag (5kg Capacity):** For carrying essentials.
- **Cane Holder:** Convenient for users needing extra support.
- **Lockable Seating:** Provides stability when seated.

INDOOR WALKER	Overall Length	Overall Width	Overall Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL	Tray SWL
RED/BLACK WAF705100RE	26"	21"	34"-38"	11"	6"	10lbs	11lbs	11lbs	11lbs
VIOLET/WHITE WAF705100VL	26"	21"	34"-38"	11"	6"	10lbs	11lbs	11lbs	11lbs

** Includes accessories, bag and tray.